











	Description	size	Ident no.	Part #
Torches made of one piece quartz glass				
	Torch with straight 8mm aerosol port		(2010069690)	800 - 05
	Torch high solids		(2010095000)	800 - 75
Semi demountable torches (Varian)				
	Torch quartz body with NS 14		(2010072900)	800 - 10
	Holding device made of PTFE		(9910057300)	800 - 48
	Injector tip made of quartz	1,8mm	(2010073000)	800 - 41
	Injector tip made of Alumina	1,8mm	(2410067590)	800 - 45
Semi demountable torch (AHF)				
	Torch quartz body			804 - 10
	Torch adaptor made of PEEK			802 - 10
	Injector tube made of Alumina	1,5mm		802 - 45
Fully demountable torch (GE)				
	Torch body & case assembly		(7910044500)	808-0805
	Outer tube		(2010095500)	808-0907
	Inner tube		(2010095300)	808-0908
	Injector made of quartz glass	1,8mm	(2010095900)	808-0906
Torch accessories				
	Bonnet (torch cap)		(2010070790)	800 - 25
	GazFit- connectors 6mm for Argon inlets (4pcs.)			GAZ-06
	GazFit- connectors 8mm for aerosol port (2pcs.)			GAZ-08
	RF coil gold-plated for Liberty II, Vista		(0110590490)	C90-1301
	RF coil gold-plated for Liberty I, Vista		(0110357990)	C90-1303
Cyclonic spray chambers made of borosilicate glass				
	without baffle with O-rings		(2010081700)	800 - 21
	with baffle and O-rings			800 - 22
	Tracey with O-rings		(2010081700)	809-0600
	Twister with O-rings		(7910043700)	809-0285
	Twister with Helix		(7910052600)	809-9199HE
Sample introduction Kits				
Varian OES radial Kit				KT-8037B
	ABC torch, fully demountable		(9910106400)	
	Twister cyclonic spray chamber with adapter			
	SeaSpray nebulizer, 2mL/min		(AR30-07-FSS2)	
AHF-Kit (complete HF-resistant)				
	Torch body			KT-AHF
	Alumina Injector	1,8mm		802 - 10
	Torch adaptor made of PEEK			802 - 45
	AHF cyclonic spray chamber made of PEEK (with baffle)			802 - 40
	PolyCon nebulizer, 2 mL/min		(AR30-07-CP2)	808 - 65
				505 - 80
Nebulizer made of borosilicate glass with UniFit sample capillary and EzyLok adapter				
	Conikal, 3mL/min	AR30-07-UC3E	(2010069500)	504 - 81
	SeaSpray, 2mL/min	AR30-07-USS2E	(2010096400)	504 - 36
	MicroMist, 0,2mL/min	AR30-07-UM02E		504 - 52
Nebulizer (HF-resistant)				
	OpalMist, 2mL/min	AR30-07-PFA2		506 - 80

Our prices can be found on our website www.ahf.de

	Description	size	Ident no.	Part #
Torches made of one piece quartz glass				
	Torch with straight 8mm aerosol port		(2010081200)	800 - 14
	Torch with aerosol port elbow		(2010090400)	800 - 16
	Torch with short outer tube for high solids		(2010094800)	800 - 13
Semi demountable torches (Varian)				
	Torch quartz body with NS 14		(2010083400)	800 - 15
	Holding device made of PTFE		(9910057300)	800 - 48
	Injector tip made of quartz	1,8mm	(2010073000)	800 - 41
	Injector tip made of Alumina	1,8mm	(2410067590)	800 - 45
Demountable Torch (AHF)				
	Torch body			804 - 20
	Holding device made of PEEK			802 - 20
	Injector made of Alumina	2,4mm		802 - 40
				802 - 46
Fully demountable torch (GE)				
	Torch body & case assembly		(7910044400)	808-0878
	Outer tube		(2010095200)	808-1231
	Inner tube		(2010095300)	808-0841
	Injector made of quartz glass	2,4mm	(2010095400)	808-0908
				808-0926
Torch accessories				
	Bonnet (torch cap)		(2010070790)	800 - 25
	GazFit- connectors 6mm for Argon inlets (4pcs.)			GAZ-06
	GazFit- connectors 8mm for aerosol port (2pcs.)			GAZ-08
	RF coil made of copper/gold for Vista AX		(0110601290)	C90-1300
	RF coil made of copper/gold for 700-ES			C90-1304
Cyclonic spray chambers made of borosilicate glass				
	without baffle with O-rings		(2010081700)	800 - 21
	with baffle and O-rings			800 - 22
	Twister mit Helix			809-0285HE
	Tracey mit O-Ringen		(2010081700)	809-0600
	Twister mit O-Ringen		(7910043700)	809-0285
Spray chamber (HF resistant)				
	made of PEEK with baffle	40ml		808 - 65
Varian OES axial Kit				
	ABC torch, fully demountable, SeaSpray nebulizer, 2mL/min			KT-8039B
	Twister cyclonic spray chamber with adapter			
AHF-Kit (complete HF resistant)				
	Torch body			KT-AHF
	Alumina Injector	1,8mm		802 - 20
	Torch adaptor made of PEEK			802 - 45
	AHF cyclonic spray chamber made of PEEK (with baffle)			802 - 40
	PolyCon nebulizer, 2 mL/min		(AR30-07-CP2)	808 - 65
				505 - 80
Nebulizers with UniFit sample capillary and EzyLok adapter				
	Conikal, 3mL/min	AR30-07-UC3E	(2010069500)	504 - 81
	Conikal, 3mL/min	AR30-1-UC3E	(201081600)	504 - 73
	SeaSpray, 2mL/min	AR30-07-USS2E	(2010096400)	504 - 36
	Slurry, 4 mL/min	AR30-07-US6E	(2010097600)	504 - 26
	Slurry, 4 mL/min	AR30-1-US6E	(2010097700)	504 - 28

Our prices can be found on our website www.ahf.de

	Description	size	Ident no.	Part #
	MS-Cones			
	Sampler made of Nickel (Ni- insert in Cu)		(G1820-65238)	900 - 71
	Skimmer made of Nickel		(G1820-65265)	900 - 72
	Sampler made of Ptlatinum(10mm Pt-insert in Cu base)		(G1820-65239)	900 - 75
	Sampler made of Platinum (18mm Pt-insert in Cu base)		(G1820-65360)	902 - 75
	Skimmer made of Platinum (7,2mm Pt-insert in Cu base)		(G1820-65237)	900 - 76
	Skimmer made of Platinum(7,2mm Pt-insert in Ni base)			902 - 76
	eSPEC Sampler made of Platinum (10mm Pt/Ir-insert)			903 - 75
eSPEC Sampler made of Platinum (18mm Pt/Ir-insert)			904 - 75	
	Torches made of quartz glass, one piece			
	Standard torch with injector tip	2,5mm	(G1820-65028)	900 - 05
	Torch with injector tip	1,5mm	(G1820-65529)	900 - 06
	Semi demountable torches			
	Torch body made of quartz glass	NS14	(G1820-65035)	900 - 08
	Injector tip made of quartz	2,0mm		900 - 25
	Injector tip made of Alumina ceramics	2,0mm		900 - 20
	Injector tip holder			900 - 30
	Accessories			
	"GazFit" connectors (4pcs.)	4/6mm		GAZ0604
	Gas connectors (2 pcs.)	4/6mm	(G1820-65027)	911 - 39
	Shield made of Platinum		(G1820-65357)	903 - 41
	Quartz sleeve for shield		(G1820-65026)	900 - 40
	Ball joint connector (quartz-cup) tubing - torch			900 - 32
	Ball joint connector (quartz-ball) chamber - tubing			900 - 36
	Ball joint adapter chamber - torch	56mm	(G1820-80237)	900 - 38
	Spray chambers made of borosilicate glass			
	Scott spray chamber with open end		(G1820-65336)	900 - 11
	Cyclonic spray chamber "Twister" with baffle	50 mL		809-0260HE
	Spray chambers made of quartz glass			
	Scott spray chamber (standard)		(G1820-65337)	900 - 13
	Scott spray chamber for Cross-Flow-Nebulizer		(G1820-65339)	900 - 14
	Spray chambers (HF resistant)			
	Spray chamber made of Teflon PFA			900 - 67
	Sample introduction kit			
	High purity kit - for semiconductor applications			Kit-4500
	Torch body made of quartz glass			902 - 10
	Injector tip made of Sapphire, ID 2,0mm			905 - 59
	Double pass spray chamber made of PFA, without KS			910 - 67
	End cap with Make-up-gas-port made of PFA	35/6		526 - 32
	MicroFlow-ST nebulizer made of PFA, 0,1mL/min			507 - 61
Nebulizers				
MicroMist nebulizer, 0,4mL/min, with UniFit			504 - 65	
Conikal nebulizer, 1mL/min, with UniFit		(G1820-65030)	504 - 70	
Cross-Flow-Nebulizer with Make-up-gas port		(G1833-65462)	501 - 93	
End cap for concentric nebulizers, PFA	35/ 6		526 - 34	

Our prices can be found on our website www.ahf.de

	Description	size	Ident no.	Part #
	MS-Cones			
	Sampler made of Nickel (Ni- insert in Copper base)		(G1820-65238)	900 - 71
	Skimmer made of Nickel (a, i, s)		(G1820-65050)	901 - 72
	Sampler (10mm Pt-insert in Copper base)		(G1820-65239)	900 - 75
	Sampler (18mm Pt-insert in Copper base)		(G1820-65360)	902 - 75
	Skimmer (7,2mm Pt-insert in Copper base)		(G1820-65237)	900 - 76
	Sampler made of Nickel (T-Mode)		(G1820-65481)	900 - 81
	Skimmer made of Nickel (T-Mode)		(G1820-65480)	900 - 82
	Skimmer 7500c made of Nickel		(G1833-65497)	902 - 82
	Skimmer 7500c made of Platinum (Nickel base)		(G1833-65092)	902 - 86
	Skimmer 7500ce made of Nickel		(G3270-65024)	903 - 82
	Skimmer 7500cs made of Platinum		(G1833-65132)	903 - 86
	Skimmer 7500cs (7,2mm Pt-insert in Nickel base)		(G3270-60106)	902 - 76
	eSPEC Sampler made of Platinum (10mm Pt/Ir)			903 - 75
	eSPEC Sampler made of Platinum (18mm Pt/Ir)			904 - 75
	eSPEC Skimmer made of Platinum, 7500cs			903 - 76
Coils				
	RF coil made of copper (silver plated)		(G1833-65432)	C90-7501
	RF coil made of copper (gold plated)			C90-7502
	Torches made of quartz glass, one piece			
	Standard torch with injector	2,5mm	(G1833-65423)	901 - 05
	Standard torch with injector	1,5mm	(G1833-65424)	901 - 06
	Standard torch with injector	1,0mm	(G1833-65425)	901 - 07
	Semi demountable torches			
	Torch quartz body			901 - 10
	Torch adaptor made of PEEK			902 - 40
	Injector made of Sapphire	2,4mm		905 - 60
	Torch accessories			
	"GazFit" connectors (4pcs.)	4/6mm		GAZ0604
	eSPEC long life "shield" plate made of Platinum	21mm	(G1833-65419)	904 - 41
	Quartz sleeve for long life "shield"	22mm	(G1833-65421)	901 - 40
	Ball joint connector chamber - torch		(G1833-65426)	901 - 38
Spray chambers made of quartz glass				
	Scott spray chamber		(G1820-65337)	900 - 13
	Scott spray chamber for Cross-Flow-Nebulizer		(G1820-65339)	900 - 14
Spray chambers (HF resistant)				
	Spray chamber made of Teflon PFA			900 - 67
	Cyclonic chamber with baffle made of PEEK	40ml		908 - 65
Sample introduction Kit				Kit-7500
	Torch body made of quartz glass			901 - 10
	Injector made of Platinum	2,0mm	(G1833-65409)	915 - 59
	Double pass spray chamber made of Teflon PFA			910 - 67
	End cap with Make-up-gas-port made of PFA	35/6		526 - 34
	MicroFlow-ST nebulizer made of PFA, 0,1mL/min			507 - 61
Nebulizer				
	Conikal, 1mL/min, with UniFit		(G1820-65030)	504 - 70
	MicroMist nebulizer, 0,4ml/min, with UniFit		(G3266-80003)	504 - 66
	Cross-flow-nebulizer with Make-up-gas port		(G1833-65462)	501 - 93

Our prices can be found on our website www.ahf.de