

	Description	Type	Ident no.	Part #
	Torches			
	Torch	Ultramass	(2010078300)	800 - 06
	Torch quartz body with NS 14	Ultramass	(2010079300)	800 - 11
	Torch	810/820	(2010100700)	800 - 34
	Injector made of quartz glass			
	ID= 1,8mm	Ultramass	(2010073000)	802 - 41
	Injector made of Alumina			
	ID= 1,8mm	Ultramass	(2410067500)	800 - 45
	ID= 2,5mm	810/820		800 - 46
	Accessories			
	Holding device made of PTFE	Ultramass	(9910057300)	800 - 48
	Holder with adjust screw for injector tip	Ultramass		800 - 47
	Bonnet (torch cap)	Ultramass	(2010070790)	800 - 25
	Sheat bypass	810/820	(2010100800)	800 - 31
	Sheat gas	810/820	(2010100900)	800 - 32
	Sample introduction Kit			
	Varian Kit , includes Torch with 1,5 mm injector Cyclonic spray chamber Twister (50 mL, with baffle), SeaSpray nebulizer (2mL/min), all adaptors			KT-8041B
	Cyclonic spray chambers made of borosilicate glass			
	without baffle		(2010081700)	800 - 21
	with baffle			800 - 22
	Twister	50 mL		809-0365HE
	Tracey, water cooled	50 mL		809-0153HE
	Spray chambers made of quartz glass			
	Scott spray chamber	810/820	(2010101000)	800 - 33
	Cyclonic spray chambers (HF resistant)			
	with baffle made of PEEK			808 - 65
	Tracey made of TFE			808 - 0963
	Nebulizers made of borosilicate glass with UniFit sample capillary and EzyLok adapter			
	Conikal, 3mL/min	AR30-1-UC3E	(2010081600)	504 - 72
	SeaSpray, 2mL/min	AR30-1-USS2E		504 - 46
	MicroMist, 0,2mL/min	AR30-1-UM02E		504 - 62
	MicroMist, 0,4mL/min	AR30-1-UM04E	(2010102100)	504 - 64
	Nebulizers (HF-resistant)			
	Polycon, 0,4mL/min	AR30-1-CP04		505 - 64
	OpalMist, 0,2mL/min	AR30-1-PFA02		506 - 62
	DuraMist, 0,4mL/min	ARG-1-DM04		505 - 24
	Microflow PFA-ST, 0,1mL/min	for oneFAST		507 - 61

Our prices can be found on our website www.ahf.de