

























	Beschreibung		Ident-Nr.	Artikel-Nr.
	Torches made of Quartz glass			
	for 6000	radial	8423 120 51141	285 - 05
	EMT for 6000	radial	8423 120 51741	285 - 13
	for 6000	DUO	8423 120 51241	285 - 12
	EMT for 6000	DUO	8423 120 51841	285 - 14
	Injector tips made of Quartz glass			
	organic	ID= 1,0mm 2 red rings	8423 120 51341	280 - 40
	aquaeuos	ID= 1,5mm 1 red ring	8423 120 51351	280 - 41
		ID= 2,0mm 1 blue ring	8423 120 51371	280 - 42
	Injector tip made of Alumina			
	HF-resistant	ID= 2,0mm	8423 120 51361	280 - 51
	Injector tips made of Quartz glass for EMT-Torch			
		ID= 1,0mm		280 - 30
		ID= 1,5mm		280 - 31
		ID= 2,0mm		280 - 36
	Injector tip made of Alumina			
	HF-resistant	ID= 2,0mm	8423 120 51961	280 - 52
	RF Coils			
	RF coil, Silver			C90-4002
	RF coil, Gold			C90-4003
	Kits for iCAP 6000 DUO (EMT)			
	Standard Kit		8423 120 51291	KT-1032A
	Organics Kit		8423 120 51201	KT-1033A
	HF-resistant Kit		8423 120 51221	KT-1035A
	Kits for iCAP 6000 Radial (EMT)			
	Standard Kit		8423 120 51191	KT-1031A
	High Solids Kit		8423 120 51131	KT-1040A
	HF-resistant Kit		8423 120 51121	KT-1039A
	Accessories			
	Periscope Window	DUO	8423 120 51591	280 - 23
	Ball joint adaptor	radial	8423 120 51151	280 - 34
	Ball joint adaptor	DUO	8423 120 51251	280 - 35
	Cyclonic spray chamber made of Borosilicate glass with Helix adapter			
	Tracey	DUO/radial	8423 120 51411	809-2665HE
	Twister	DUO/radial	8423 120 51311	809-2951HE
	Spray chamber (HF-resistant)			
	Tracey TFE	radial	8423 120 51451	809-9198
	Tracey TFE	DUO	8423 120 51461	809-0963
	Cyclonic, PEEK	DUO		278 - 23
	Nebulizers made of Borosilicate glass			
	Conikal, 2mL/min, with UniFit	AR30-07-UC2E	8423 120 51431	504 - 80
	Aerosalt		8423 120 51331	500 - 65
	Nebulizers (HF resistant)			
	VeeSpray, Ceramic	AR30-07-CV6	8423 120 51321	516 - 30
	OpalMist, 2mL/min	AR30-07-PFA2	8423 120 51471	506 - 80
	Burgener Mira Mist PTFE		8423 120 51471	508 - 54

Our prices can be found on our website www.ahf.de

	Description	Ident-no.	Part-no..
Semi demountable torches in ceramic-bodies (incl. injectors)			
	Torch for TRACE, AtomScan, IRIS AP	(12643203)	275 - 10
	Dual-view-model for TRACE DV (incl. sleeve)	(13856200)	275 - 12
Semi demountable torches in ceramic-bodies (without injectors)			
	Torch for TRACE, AtomScan, IRIS AP	(12644003)	275 - 101
	Dual-view-model for TRACE DV	(12644005)	275 - 121
Fully demountable torches (ABC-Torch)			
		DUO	808-0692
	Quartz outer tube		808-0688
	Quartz inner tube, high flow		808-0972
	Tapered quartz injector tube, ID= 1,8mm		808-0690
		Axial	808-0704
	Quartz outer tube		808-0707
	Quartz inner tube, high flow		808-0972
	Tapered quartz injector tube, ID= 1,8mm		808-0690
Injector tubes made of quartz			
	ID= 1,0mm 2 red rings	for organic samples	(12540701) 250 - 40
	ID= 1,5mm 1 red ring	for aqueous samples	(12540700) 250 - 41
	ID= 2,0mm 1 blue ring		(13650201) 250 - 42
	ID= 2,4mm 2 blue rings		(13650202) 250 - 43
Injectors made of Alumina			
	ID= 1,5mm		(13643601) 257 - 41
	ID= 1,8mm		(13643602) 257 - 42
Injectors made of Sapphirre			
	ID= 1,5mm		254 - 41
	ID= 2,0mm		254 - 42
	ID= 2,5mm		254 - 43
Accessories			
	Anti arc cap		(08520700) 275 - 25
	Small bonnet (flat)		(08520701) 275 - 26
	DUO torch sleeve made of quartz		(13856300) 250 - 23
	Torch adapter for cyclonic spray chamber		(14174300) 275 - 20
Spray chamber made of borosilicate-glass			
	Cyclonic for ICAP 61E trace analyser (siphon drain)		(13462000) 275 - 15
	Cyclonic for IRIS and Atomscan series (siphon drain)		(13612000) 275 - 16
	Cyclonic for IRIS and Atomscan series (pump drain)		(13966200) 275 - 17
Caclonic spray chambers made of PEEK (HF resistant)			
	vertical circulation		278 - 23
	horizontal circulation		278 - 63
Nebulizer made of borosilicate glass (with UniFit and EzyLok connectors)			
	Conikal, 1mL/min	AR30-07-UC1E	(1210000) 504 - 79
	Conikal, 2mL/min	AR30-07-UC2E	(13918400) 504 - 80
	SeaSpray, 2mL/min	AR30-07-USSE	(14186300) 504 - 36
Nebulizer (HF-resistant)			
	Cross-flow-nebulizer (HF-resistant) with Pt-capillary and sapphirre nozzle		
	in PTFE-endcap for spray chamber ID= 35mm		(12698900) 500 - 92
	Polycon made of Polyimide, 2mL/min		(14189900) 505 - 80
	Burgener Mira Mist made of PTFE, 0,2 - 2,5mL/min		(13638800) 508 - 54

Our prices can be found on our website www.ahf.de

	Description	Ident-no.	Part-no.
Torches			
	Torch, injectortube 1,0 mm for TJA 975, 9000, 61&1100	(19040054)	250 - 05
	Torch for organic samples, inj. 1,0 mm for IRIS RP	(13913700)	252 - 08
	Torch bonnet	(19040075)	250 - 25
Semi demountable torches in ceramic bodies (incl. injectors)			
	High/low-flow-model for IRIS RP, PolyScan	(12643202)	275 - 06
	High-flow-model for IRIS RP, PolyScan	(12643201)	275 - 07
	High-flow-model for IRIS RP, PolyScan (for organics)	(12512401)	275 - 08
Semi demountable torches in ceramic bodies (without injectors)			
	High/low-flow-model for IRIS RP, PolyScan	(12644002)	275 - 061
	High-flow-model for IRIS RP, PolyScan	(12644001)	275 - 071
Fully demountable torches (ABC-Torch)			808-0696
	Quartz outer tube		808-0691
	Quartz inner tube, high flow		808-0972
Injector tubes made of quartz			
	ID= 1,0mm	2 red rings for organic samples	(12540701) 250 - 40
	ID= 1,5mm	1 red ring for aqueous samples	(12540700) 250 - 41
	ID= 2,0mm	1 blue ring standard	(13650201) 250 - 42
	ID= 2,4mm	2 blue rings for better DL's	(13650202) 250 - 43
Injectors made of Alumina ceramics			
	ID= 1,5mm		257 - 41
	ID= 1,8mm		257 - 42
Injectors made of Sapphire			
	ID= 1,5mm		254 - 41
	ID= 2,0mm		254 - 42
	ID= 2,5mm		254 - 43
Accessories			
	Anti arc cap	(08520700)	275 - 25
	Small bonnet (flat)	(08520701)	275 - 26
Spray chambers made of borosilicate-glass			
	Cyclonic for IRIS Advantage Floor Model (pump drain)	(13966300)	275 - 14
	Cyclonic for IRIS RP-spectrometer (pump drain)	(14054200)	275 - 18
	Cyclonic for IRIS RP-spectrometer (siphon drain)	(13610500)	275 - 19
	Impact bead spray-chamber	(13921900)	275 - 21
Spray chamber (HF resistant)			
	Cyclonic spray chamber made of PEEK with baffle		278 - 65
Nebulizer made of borosilicate glass (with EzyFit and EzyLok side arm)			
	Conikal, 1mL/min	AR30-07-FC1E	(1210000) 503 - 79
	Conikal, 2mL/min	AR30-07-FC2E	(13918400) 503 - 80
	SeaSpray, 2mL/min	AR30-07-FSSE	(14186300) 503 - 36
Nebulizer (HF-resistant)			
	Cross-Flow-Nebulizer (HF-resistant) with Pt-capillary and saphirre nozzle in PTFE-endcap for spray chamber ID= 35mm		(12698900) 500 - 92
	Polycon made of Polyimide, 2mL/min		(14189900) 505 - 80
	Burgener MiraMist made of PTFE, 0,2 - 2,5mL/min		(13638800) 508 - 54

Our prices can be found on our website www.ahf.de