

PFA CONTAINER ultra pure

Price List 2010

high purity polymer
 no catalysts or stabilizers traces
 no blanks and memory effects


extremely smooth and hydrophobic surface
 no adsorption and desorption
 non-wetting, non-sticking

easy and quick cleaning
 superior thermal stability -200° up to +260°C
 sterilization with autoclave up to 200°C




D [mm]	H [mm]	volume	part no.	SP	€/ set
--------	--------	--------	----------	----	--------

Sample cups		inside conical bottom handling with tweezers possible with calibration marks in mL extra transparent wall shows content				
13,6	24	fits in autosampler rack	1,5 mL	T43 - 002	2	31,50
14,0	24		1,5 mL	T43 - 022	2	29,50
13,6	36	fits in autosampler rack	2,5 mL	T43 - 003	2	33,80
14,0	36		2,5 mL	T43 - 023	2	32,00
14,0	52		4,0 mL	T43 - 024	2	46,50



Caps, PFA						
18	5	tight sealing function		T43 - 918	2	29,60
16	10	for dust protection		T43 - 916	2	39,50



Sample vials, without closure, standing for ICP-Autosampler						
15	50	flat bottom	4 mL	T33 - 505	1	8,70
15	43	round bottom	6 mL	T33 - 506	1	11,20
22	55	round bottom	12 mL	T33 - 515	1	28,90
27	59	flat bottom	20 mL	T33 - 522	1	11,90
27	59	round bottom	25 mL	T33 - 525	1	18,60
		snap-on cap, PFA	for 4mL	T33 - 500	2	46,50
			for 12mL	T33 - 510	2	29,60
			for 20mL	T33 - 520	2	52,50



Sample vials, with snap-on caps, standing for ICP-Autosampler						
15	50	flat bottom	4 mL	T33 - 555	1	13,35
27	59	flat bottom	20 mL	T33 - 560	1	17,15



Sample vials, with threaded closure, standing for ICP-Autosampler						
22	55	flat bottom	12 mL	T33 - 512	1	12,30

Wide neck volumetric flasks with screw cap (for ESI SC-Autosampler)						
147	neck ID= 16mm		50 mL	T87 - 050	1	96,30
177	neck ID= 16mm		100 mL	T87 - 100	1	124,00



SP 1 = per piece; SP 2 = per 10 pieces

all prices ex works (Tübingen)