

high purity polymer
no catalysts or stabilizers traces
no blanks and memory effects


extremely smooth and hydrophobic surface
no adsorption and desorption
non-wetting, non-sticking

easy and quick cleaning
superior thermal stability -200° up to +260°C
sterilization with autoclave up to 200°C




| D [mm] | H [mm] | volume | part no. | SP | € / set |
|--------|--------|--------|----------|----|---------|
|--------|--------|--------|----------|----|---------|

| Sample cups | | | | | |
|--|----|--------------------------|---------------|-----------|----------------|
| inside conical bottom handling with tweezers possible with calibration marks in mL extra transparent wall shows content | | | | | |
| 13,6 | 24 | fits in autosampler rack | 1,5 mL | T43 - 002 | 2 31,50 |
| 14,0 | 24 | | 1,5 mL | T43 - 022 | 2 29,50 |
| 13,6 | 36 | fits in autosampler rack | 2,5 mL | T43 - 003 | 2 33,80 |
| 14,0 | 36 | | 2,5 mL | T43 - 023 | 2 32,00 |
| 14,0 | 52 | | 4,0 mL | T43 - 024 | 2 46,50 |



| Caps, PFA | | | | | |
|------------------|----|------------------------|--|-----------|----------------|
| 18 | 5 | tight sealing function | | T43 - 918 | 2 28,40 |
| 16 | 10 | for dust protection | | T43 - 916 | 2 41,70 |



| Sample vials, without closure, standing for ICP-Autosampler | | | | | |
|--|----|--------------|--------------|-----------|----------------|
| 15 | 50 | flat bottom | 4 mL | T33 - 505 | 1 8,70 |
| 15 | 43 | round bottom | 6 mL | T33 - 506 | 1 11,20 |
| 22 | 55 | round bottom | 12 mL | T33 - 515 | 1 28,90 |
| 27 | 59 | flat bottom | 20 mL | T33 - 522 | 1 11,90 |
| 27 | 59 | round bottom | 25 mL | T33 - 525 | 1 18,60 |
| snap-on cap, PFA | | | for 4mL | T33 - 500 | 2 46,50 |
| | | | for 12mL | T33 - 510 | 2 29,60 |
| | | | for 20mL | T33 - 520 | 2 52,50 |



| Sample vials, with snap-on caps, standing for ICP-Autosampler | | | | | |
|--|----|-------------|--------------|-----------|----------------|
| 15 | 50 | flat bottom | 4 mL | T33 - 555 | 1 13,35 |
| 27 | 59 | flat bottom | 20 mL | T33 - 560 | 1 17,15 |



| Sample vials, with threaded closure, standing for ICP-Autosampler | | | | | |
|--|----|-------------|--------------|-----------|----------------|
| 22 | 55 | flat bottom | 12 mL | T33 - 512 | 1 11,40 |

| Wide neck volumetric flasks with screw cap (for ESI SC-Autosampler) | | | | | |
|--|---------------|--|---------------|-----------|-----------------|
| 147 | neck ID= 16mm | | 50 mL | T87 - 050 | 1 96,30 |
| 177 | neck ID= 16mm | | 100 mL | T87 - 100 | 1 124,00 |



SP 1 = per piece; SP 2 = per 10 pieces

all prices ex works (Tübingen)

