


## DESICCATORS

Dry and clean storage of sensitive materials and objects



### PROTECT WHATS VALUABLE

Desiccators enable you to safely store and transport air- or moisture-sensitive materials or cleanroom-manufactured components that have to remain dust-free. The desiccators are convenient in handling and compact in design to save lab space. Various sizes and designs are available to fit your requirements.

### PROTECTION AGAINST ENVIRONMENTAL INFLUENCE

- :: Dust
- :: Humidity
- :: Oxygen



## AUTO DESICCATORS

**Perfect for long-term storage of sensitive samples**

- :: Stable and continuous dehumidification by polymer electrolyte membrane
- :: Keeps humidity at constant level
- :: All models equipped with digital hygrometer



## STANDARD DESICCATORS

**For storage of sensitive materials**

- :: Stacking possible
- :: Includes shelves and pan for silica gel
- :: Smoked glass type desiccators available for preventing deterioration or discoloring
- :: Optional with one-touch-door (open and close with one hand)



## MINI DESICCATORS

**For safe transport or long-term storage at low humidity**

- :: Stable and uniform quality by molded design
- :: Includes shelves and pan for silica gel
- :: Optional with handle
- :: Space-saving shelving



## VACUUM DESICCATORS

**For storage in inert atmosphere, packaging tests, degassing, etc.**

- :: Air-tightness through high-lever bonding
- :: Made of heavy-walled acrylic polymer
- :: Super transparent
- :: 2 valves for vacuum/gas connection
- :: Pressure gauge included

### **Glove boxes:**

For working in inert gas surrounding. Please ask our team.